Argonne

www.anl.gov

ENVIRONMENTAL SCIENCE • PLANETARY SCIENCE

CHEMICAL SCIENCE •

ATOMIC PHYSICS • MOLECULAR PHYSICS • X-RAY INSTRUMENTATION • MATERIALS SCIENCE • BIOLOGICAL SCIENCE •

Advanced Photon Source Bldg. 401/Rm A4115 Argonne National Laboratory 9700 S. Cass Ave Argonne, IL 60439 USA apsinfo@aps.anl.gov • www.aps.anl.gov

ADVANCED PHOTON SOURCEA New Tool for Structure Studies at High Pressures with a DAC

A scanning angle energy-dispersive x-ray diffraction (SA-EDXD) technique employing an energy-dispersive solid-state detector and white synchrotron radiation has been developed for high-pressure structure studies using a diamond anvil cell (DAC). The main feature of the technique is well-defined collimation of the diffracted beam path into the detector, which improves the signal-to-noise ratio significantly compared to routine monochromatic angle-dispersive powder diffraction. This is particularly useful and essential for materials with low-scattering cross section and/or amorphous and liquid diffraction/scattering studies using a DAC.

The diamond anvil cell is one of the most common devices used for high-pressure structure studies. Due to the thin sample used in a

DAC at high pressure, the background from the diamond and from pressure-transmitting media is often significant for routine angle-dispersive powder diffraction. Combining a solid-state point detector with a scanning diffraction angle, one can establish good collimation to measure only the diffraction signals from a sample (Fig. 1). A group at the High Pressure Collaborative Access Team (sector 16 at the APS) has extended the technique—making it suitable for DAC experiments—by employing high-precision stages with a sphere-of-confusion that is less than 5 μ m. Together with software developed for automation control and data analysis, this technique is now available for APS general users. The diffraction intensity is sorted to 4,000 channels by the energy-dispersive detector. When the multichannel energy data is converted to constant-energy angle-dispersive data, only a coarse 20 scan is required to achieve the regular angle-

dispersive diffraction resolution. One can achieve an angle-dispersive diffraction pattern at certain energies (and energy-dispersive diffraction pattern at other angles) for structure analysis and refinement.

The SA-EDXD technique provides an important new tool for amorphous and liquid structure studies under extreme conditions. The well-defined collimation ensures the necessary signal-to-noise ratio for weak signals from amorphous and liquid materials in DAC. The use of white beam up to 110 keV covers a large Q-range (>20 Å⁻¹). The intensity data collected by the SA-EDXD technique are a function of energy and angle, from which one can obtain angle-dispersive data from different energies, and energy-dispersive data at different angles as well. The structure factor, S(Q), can be calculated from the angle-dispersive pattern at different energies (Fig. 2), while each one has different Q-range. A combined S(Q) can be obtained for a large Q-range. The work for partial structural factors is in progress.

Contact: Wenge Yang (wyang@hpcat.aps.anl.gov)

See: Y. Wang et al., J. App. Cryst. (2004) **37**, 947-956.

This work is supported by the U.S. DOE-Basic Energy Sciences, DOE-National Nuclear Security Administration, NSF, DOD-TACOM, and the W.M. Keck Foundation. Use of the Advanced Photon Source is supported by the U.S. Department of Energy, Office of Science, Office of Basic Energy Sciences, under Contract No. No. DE-AC02-06CH11357.




Fig. 1. Schematic of data acquisition with the SA-EDXD technique.




Fig. 2. Nine angle-dispersive data (a), with energies from 24 keV to 74 keV, are used to obtain the full structural factor S(Q).

Argonne is a U.S. Department of Energy laboratory managed by UChicago Argonne, LLC

CALL FOR PROPOSALS

At the Advanced Photon Source, our door is open to experimenters from all scientific disciplines whose research requires the highest brilliance hard x-ray beams in the Western Hemisphere.

General-user proposals for beam time during Run 2007-3 are due by July 13, 2007.

Information on access to beam time at the APS is at http://www.aps.anl.gov/user/beamtime/get_beam.html or contact Dr. Dennis Mills, DMM@aps.anl.gov, 630/252-5680.

A U.S. Department of Energy laboratory managed by UChicago Argonne, LLC The Advanced Photon Source is funded by the U.S. Department of Energy, Office of Science, Office of Basic Energy Sciences

lightsources.org

www.lightsources.org



