## Microsymposium

## MS54.003

## Electron microscopy and functional analysis of recombinant innexin gap junctions

<u>A. Oshima<sup>1</sup>, T. Matsuzawa<sup>2</sup>, K. Murata<sup>3</sup>, K. Nishikawa<sup>1</sup>, Y. Fujiyoshi<sup>1</sup></u>

<sup>1</sup>Nagoya University, Cellular and Structural Physiology Institute, Nagoya, Japan, <sup>2</sup>Kyoto University, Department of Biophysics, Graduate School of Science, Kyoto, Japan, <sup>3</sup>National Institute for Physiological Sciences, Okazaki, Japan

Innexin is a molecular component of invertebrate gap junctions, which have an important role in neural and muscular electrical activity in invertebrates. Although the structure of vertebrate connexin26 was revealed by X-ray crystallography [1], the structure of innexin channels remains poorly understood. To study the structure of innexin gap junction channels, we expressed and purified Caenorhabditis elegans innexin-6 (INX-6) gap junction channels, and characterized their molecular dimensions and channel permeability using electron microscopy (EM) and a fluorescent dye transfer assay, respectively [2]. Negative-staining and thin-section EM of isolated INX-6 gap junction plaques revealed a loosely packed hexagonal lattice. We performed single particle analysis of purified INX-6 channels with negative-staining and cryo EM. Based on the negative-stain EM images, the class average of the junction form had a longitudinal height of 220 Å, a channel diameter of 110 Å in the absence of detergent micelles, and an extracellular gap space of 60 Å, whereas the class average of the hemichannels had diameters of up to 140 Å in the presence of detergent micelles. Cryo EM images revealed rotational peaks that could be related to the INX-6 subunits. Structural analysis of the reconstituted INX-6 channels with single particle analysis and electron tomography suggested that the oligomeric number of the INX-6 channel was distinct from that of the dodecameric connexin channel. Dye transfer experiments indicated that the INX-6 channels have a characteristic oligomer component that differs from that in connexin gap junction channels.

[1] S. Maeda, S. Nakagawa et al., Nature, 2009, 458, 597–602, [2] A. Oshima, T. Matsuzawa et al., J. Biol. Chem., 2013, 288, 10513-10521

Keywords: Gap junctions, Innexin, Electron microscopy