## **Keynote Lectures**

M. Carcelli, *CrystEngComm*, **2005**, *7*, 527–537. [6] A. Bacchi, G. Cantoni, M.R. Chierotti, A. Girlando, R. Gobetto, G. Lapadula, P. Pelagatti, A. Sironi, M. Zecchini, *CrystEngComm*, **2011**, DOI:10.139/c0ce00816h

Keywords: supramolecular, inclusion, organometallic

## **KN13**

Acta Cryst. (2011) A67, C10

Bacterial pathogenesis and peptidoglycan degradation machines Juan A. Hermoso, Department of Crystallography and Structural Biology, Instituto de Química Física "Rocasolano", CSIC, Madrid. (Spain). E-mail: xjuan@iqfr.csic.es

The bacterial cell wall is comprised of cross-linked strands of peptidoglycan (PG), which encase the entire cytoplasm. A healthy cell wall is critical for survival of bacteria and serves as a docking station for bacterial surface proteins, some of them representing key players in adhesion, colonisation and virulence. Therefore cell wall remodeling is critical in host-pathogen interactions, cell division, virulence, PG recycling and antibiotics resistance. Choline-binding proteins (Cbp) are pneumococcal surface proteins that have received considerable attention because of their versatility, and their sophisticated role in the interaction with host proteins. The three-dimensional structure of Pce (70 KDa), in complex with the reaction product and choline analogs has been solved [1]. We have showed that Pce hydrolyses PAF, a potent lipidic first messenger of inflammatory processes. Besides, the structural analysis indicated that Pce selectively remodels the bacterial surface impairing the ability of host proteins to efficiently bind the bacteria, and would provide a mechanism for pneumococci escaping attack by the host defense system [1].

CbpF (38 KDa) has been solved [2] demonstrating that CbpF inhibit the activity of autolysin LytC and providing the first example of a regulatory system to tune the activity of an autolysin. LytC is involved in the virulence mechanism of fratricide. Pneumococci that are competent for natural genetic transformation kill and lyse noncompetent sister cells or members of related species that are present in the same environment. This phenomenon has been termed fratricide. We have reported the functional characterization of the key effector of pneumococcal fratricide CbpD [3] and the crystal structure of LytC in a ternary complex with choline and a PG fragment [4] that explain the activation of LytC by CbpD in fratricide and provide the first structural insights into the critical and central function that LytC plays in pneumococcal virulence.

During homeostasis, including growth, cell wall is simultaneously biosynthesized and degraded. Lytic transglycosylases (LTs) initiate the degradative events on cell wall. The products of LTs are internalized to the cytoplasm, and hydrolised by glucosaminidase NagZ and amidase AmpD. The reaction products of AmpD play roles in both PG recycling events and in an induction event that leads to the expresion of  $\beta$ -lactamase, a key  $\beta$ -lactam antibiotic resistance enzyme. Despite the 3D structure for AmpD enzyme was known by NMR techniques, we have performed crystal structure determination of some LTs and AmpD. Unexpectedly we have observed that AmpD undergoes an activation mechanism from an inactive form (that determined by NMR) to an active form we solved. Changes produced in the activation process are among the largest structural rearrangements ever reported for a single domain protein.

[1] J. A. Hermoso, L. Lagartera, A. González, M. Stelter, P. García, M. Martínez-Ripoll, M., J. L. García, M. Menéndez, *Nature Structural & Molecular Biology* **2005**, *12*, 533-538. [2] R. Molina, A. González, M. Stelter, I. Pérez-Dorado, R. Kahn, M. Morales, S. Campuzano, N. E. Campillo, S. Mobashery, S., J. L. García, P. García and J. A. Hermoso *EMBO reports* **2009**, *10*, 246-

251. [3] V. Eldholm, O. Johnsborg, D. Straume, H. S. Ohnstad, K. H. Berg, J. A. Hermoso and L. S. Håvarstein *Molecular Microbiology* **2010**, *76(4)*, 905-917. [4] I. Pérez-Dorado, A. González, M. Morales, R. Sanles, W. Striker, W. Vollmer, S. Mobashery, J. L. García, M. Martínez-Ripoll, P. García and J. A. Hermoso *Nature Structural & Molecular Biology* **2010**, *17*, 576-582.

Keywords: Biocrystallography, Disease, Pathogenesis

## **KN14**

Acta Cryst. (2011) A67, C10

Coherence in crystallography for imaging materials and biology

Ian Robinson, a,b aLondon Centre for Nanotechnology, University College, London (UK). bResearch Complex at Harwell, Harwell-Oxford Campus, Didcot (UK). E-mail: i. robinson@ucl.ac.uk

In crystallography, we consider nanometre-sized crystals to be a new frontier of opportunity to tailor physical properties using 'size' as a control variable. This general view of crystallography "beyond the unit cell" opens a wide field of opportunity for methods development. However, when we think about nanostructures, we must reconsider the standard bulk concepts of lattices and crystal defects. Changes here provide nanomaterials with new and exciting properties. This lecture is about coherence-based methods for the determination of crystal structures using coherence. One example is micron-sized ZnO crystals, attached by bonding to a SiO<sub>2</sub> substrate, which showed internal strain arising from accidental damage during manipulation. Use of more than one Bragg peak from the same crystal allows components of the full strain tensor to be mapped inside the crystal.

These new crystallographic methods have a fundamental need for beam coherence, so benefit directly from 3rd generation synchrotron sources. The coherence leads to interference fringes in the diffraction patterns of sufficiently small crystals. When the fringes are measured using a fine-pixel detector the data can be oversampled beyond their Nyquist/Shannon frequency. To invert the diffraction, we then solve the crystallographic 'phase problem' using a support-constrained HIO algorithm. This leads to quantitative three-dimensional maps of the density of the crystal with a real-space phase, which is interpreted as the deformation of a crystal from its equilibrium lattice spacing.

We have also used the methods to examine twin domain structures within crystals as well as real-space phase domains due to ordering. In biology, we have investigated collagenous tissue, which shows an analogous phase domain structure. For extended objects such as these, we are applying the principles of 'ptychography', in which the crystallographic phase information emerges from overlaps between coherent diffraction patterns.

Keywords: phasing, nanomaterials, ptychography

## **KN15**

Acta Cryst. (2011) A67, C10-C11

Structure-Property Relationships of Inclusion Compounds
Len Barbour, Department of Chemistry, University of Stellenbosch,
Stellenbosch (South African). E-mail: ljb@sun.ac.za

Multicomponent crystals based on the encapsulation of guest molecules by suitable hosts may exhibit a number of interesting properties and structural phenomena. For example, guest-templated