Exploring the dynamic of PSII at room temperature by simultaneous femtosecond X-ray spectroscopy and dffraction

Photosystem II (PSII) is a membrane protein that catalyzes photo oxidation of water into dioxygen, one of nature's most fascinating and important reactions. The water-splitting reaction takes place at the oxygen evolving complex (OEC), through five intermediate S-states (S_0 to S_4)¹, where S_1 is the dark-stable state and S_3 is the last semi-stable state before O-O bond formation and O_2 evolution. The structure of PSII in the dark state has been solved by X-ray diffraction and X-ray free electron laser (XFEL) providing information of the structural geometry of the Mn_4CaO_5 cluster in OEC. However, to fully understand the O-O bond formation mechanism, elucidating the structures of the OEC in the different S-states is essential.

In our recently published study, we report high-resolution structures of PSII at room temperature using XFEL coupling with X-ray spectroscopy under different illumination conditions. This enables us to gain new insights about the dynamic changes in the structure of the Mn_4CaO_5 cluster as well as the ligands and the bound water molecules. We further use ammonia as water analogue to investigate the water-binding site(s) and discriminate between mechanisms proposed in literature. We discuss the precise role of the water bound to OEC in electron transfer and the water-splitting reaction.

LouiseLassalle Iris Young Mohamed Ibrahim Ruchira Chatterjee Sheraz Gul Franklin Fuller Aaron Brewster	USA USA Germany USA USA USA USA	Lawrence Berkeley National lab Lawrence Berkeley National lab Humboldt University Lawrence Berkeley National lab Lawrence Berkeley National lab Lawrence Berkeley National lab
Nicholas Sauter	USA	Lawrence Berkeley National lab
Athina Zouni	Germany	Humboldt University
Jan Kern	USA	Lawrence Berkeley National lab
Vittal Yachandra	USA	Lawrence Berkeley National lab
Junko Yano	USA	Lawrence Berkeley National lab